

Absorbent

Paper Points

Efficient Moisture Absorption,
Precise Oral Health Protection



RAPID DRYING

TOUGH MATERIAL

STRICT STERILITY

CLEAR IDENTIFICATION

www.tiadent.com



In the intricate world of dental treatments, every step is crucial to the success of the procedure. Absorbent Paper Points, the trusted ally of dentists, safeguard your oral health with outstanding performance.

Dimensions:

The product dimensions comply with the corresponding charts in YY 0711-2009.


High Absorbency:

The paper point exhibits excellent absorbency, with the liquid absorbed by the paper point reaching at least 10mm above the liquid surface.

Degradation Resistance:

When the paper point is immersed in water at 40°C ± 1°C for 10 minutes and then removed with tweezers, it does not disintegrate.

Taper 02	015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055 060 ,070 ,080 ,015-040 , 045-080
Tornado Blue Taper 04	015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055 060 ,070 ,080 ,015-040 , 045-080
Tornado Blue Taper 06	015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055 060 ,070 ,080 ,015-040 , 045-080
Tornado Gold	TF1 ,TF2 ,TF3 ,TF4 ,TF5 ,TF1-TF3 ,TF4-TF5
Viera Gold	TV20 ,TV25 ,TV35 ,TV45 ,TV25-TV45
Premium Blue	T25 ,T40 ,T50 ,T25-T50
Premium Gold	TX1 ,TX2 ,TX3 ,TX4 ,TX5 ,TX1-TX3 ,TX4-TX5



TiaDent Inc.
1627 Cortex Midtown Pl Rosenberg, TX 77471

+ 13184800797