

Gutta Percha

Points

Precise and reliable root canal filling solution



ACCURATE SIZE

CLEAR X-RAY IMAGING

PURE INGREDIENTS

CLEAR IDENTIFICATION

www.tiadent.com



Accurate Size

Our Gutta Percha points are produced with advanced technology and strict audit standards, ensuring dimensions are accurate to $\pm 0.01\text{mm}$ for a reliable fit.



Clear X-Ray Imaging

The Gutta Percha points are designed to provide clear visibility under X-rays, allowing for easy monitoring of the root canal filling.



Clear Identification

Each point is clearly labeled with standard specifications to ensure easy identification during procedures.



Pure Ingredients

Manufactured using high-purity raw materials, our Gutta Percha points provide enhanced therapeutic effects for optimal clinical results.

Taper 02 015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055
060 ,070 ,080 , 015-040 , 045-080

Tornado Blue Taper 04 015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055
060 ,070 ,080 ,015-040 , 045-080

Tornado Blue Taper 06 015 ,020 ,025 ,030 ,035 ,040 ,045 ,050 , 055
060 ,070 ,080 ,015-040 , 045-080


Tornado Gold TF1 ,TF2 ,TF3 ,TF4 ,TF5 ,TF1-TF3 ,TF4-TF5

Viera Gold TV20 ,TV25 ,TV35 ,TV45 ,TV25-TV45

Premium Blue T25 ,T40 ,T50 ,T25-T50

Premium Gold TX1 ,TX2 ,TX3 ,TX4 ,TX5 ,TX1-TX3 ,TX4-TX5

LENGTH MARKED



TiaDent Inc.

1627 Cortex Midtown Pl Rosenberg, TX 77471

+ 13184800797