

TIA DIAMOND

BURS

Uniform Diamond Coating

Dynamic Raw Material Selection

Optimal Embedding Of Diamond Particles

Excellent Concentricity And Durability

● Standard
● Fine
● Coarse
● Ultra-fine



www.tiadent.com

GRAPHIC TIA DIAMOND BURS

173/014 — ISO No. (Shapes/Max. Diameter)

10.0/24 — Working Length / Full Length

TIA-TF-11 — Order NO.

The color mark for showing the thickness of burs

Apical diameter (1/10mm)	Rotary speed (rpm)	Maximum speed (rpm)
*08-010	75,000-150,000	450,000
012-014	60,000-110,000	450,000
016-018	45,000-88,000	450,000
021-023	40,000-75,000	300,000
025-027	30,000-65,000	160,000
029-031	25,000-56,000	140,000
033-040	22,000-45,000	120,000
042-050	20,000-37,000	95,000
055-060	17,000-32,000	80,000
065-080	13,000-26,000	60,000


ADVANTAGE

- Ⓧ Easy to place
- Ⓧ Easy to retreat, like gutta percha
- Ⓧ Easy to remove
- Ⓧ Dense 3D obturation
- Ⓧ Easy to create post space
- Ⓧ Latex-free gutta percha

MODEL LIST



VENEER	Model	ISO No.
	TIA-EX-51	522/029
	TIA-EX-52	552/016
	TIA-EX-53	552/018
	TIA-EX-54	552/021



TiaDent Inc.
1627 Cortex Midtown Pl Rosenberg, TX 77471

+ 13184800797