

TIA TORNADO GOLD PATHGLIDER

**LOW PROBABILITY OF
LEDGE FORMATION**

**LOW RISK OF APICAL FORAMEN
DISPLACEMENT**

**EFFECTIVELY NAVIGATE TIGHT
CANAL CURVATURES**



PRIME

Size	Taper	Length	Torque	Speed	Packing
● TTG13	02	21,25,31mm	2-3.0N/cm	300rpm	6pcs/box
○ TTG16	02	21,25,31mm	2-3.0N/cm	300rpm	6pcs/box
● TTG19	02	21,25,31mm	2-3.0N/cm	300rpm	6pcs/box

Cross-Section

