


TIA HAND

FILES

Endodontic Instruments

ISO standard, suitable for many types of root canals; For removal of pulp tissue and debris, as well as preliminary root canal preparation.




www.tiadent.com




TIA K FILES

Material : Stainless Steel/NiTi
Length : 21,25,28,31 mm




TIA H FILES

Material : Stainless Steel/NiTi
Length : 21,25,28,31 mm



TIA REAMERS

Material : Stainless Steel/NiTi
Length : 21,25,28,31 mm



TIA SPREADERS

Material : Stainless Steel/NiTi
Length : 21,25,28,31 mm



TIA PLUGGERS

Material : Stainless Steel/NiTi
Length : 21,25,28,31 mm



TIA C FILES

Material : Stainless Steel
Length : 21,25,31 mm



TIA NERVE BROACHES

Material : Stainless Steel
Length : 21,25,28 mm



DESCRIPTION:

Thread hole to prevent dropping

A thread allows a chain or string to pass through it to prevent the reamer from dropping when during use.

Icons and numbers allow easy identification

The top of the grip is marked with an easy-to-identify icon with a number.

Anti-slip molded rings

The grip is molded to have anti-slip rings

Colors allow easy identification


The grips are grouped into different colors according to their thicknesses

Accurate dimensional machining

.02 taper meets ISO standards

Curved design for better holding

Easy-to-hold shape that fits the fingers



TiaDent Inc.

1627 Cortex Midtown Pl Rosenberg, TX 77471

+ 13184800797